

## Experiment No. 1

### Estimation of sodium carbonate and sodium hydrogen carbonate present in a mixture.

Date:

Principle:

Neutralisation of  $\text{Na}_2\text{CO}_3$  solution by strong acid (HCl) occurs in two steps:



The acid base indicator is to be so selected that its  $\text{P}^{\text{H}}$  range for the colour change coincides with the sudden sharp change of  $\text{P}^{\text{H}}$  at the equivalence point. So at the first neutralisation point ( $\text{P}^{\text{H}} = 8.3$ ), phenolphthalein ( $\text{P}^{\text{H}}$  range 3.1-10) shows its colour change from pink to colourless. At this stage  $\text{Na}_2\text{CO}_3$  consumes only half the amount of HCl required for complete neutralisation. If methyl orange ( $\text{P}^{\text{H}}$  range = 3.1 to 4.4) is added to this titrated solution and the titration with HCl is continued upto second equivalence point, then this titre value corresponds to the amount of HCl required to convert  $\text{NaHCO}_3$  to NaCl (i.e.  $\text{NaHCO}_3$  derived from  $\text{Na}_2\text{CO}_3$  plus the amount of  $\text{NaHCO}_3$  present in the original mixture).

1000ml 1(N) HCl  $\equiv$  53 g of  $\text{Na}_2\text{CO}_3 \equiv$  84 g of  $\text{NaHCO}_3$

#### Experimental Data

**Table-1:** Standardization of supplied NaOH solution with the supplied standard (N/20) oxalic acid solution

No of titration	Volume of oxalic acid solution (ml)	Volume of NaOH solution required (ml)	Mean volume of NaOH solution (ml)	Strength of NaOH solution
1	25	22.1	22.1	1.1312 (N/20)
2	25	22.1		

#### Calculation:

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$\text{Strength of NaOH solution (S}_2) = \frac{25 \times \text{N}/20}{22.1} = 1.1312 \text{ (N/20)}$$

$V_1$  = Vol. of oxalic acid solution

$S_1$  = Strength of oxalic acid solution

$V_2$  = Vol. of NaOH solution

**Table-2:** Standardization of supplied HCl solution with the standardized (N/20) NaOH solution

No of titration	Volume of HCl solution (ml)	Volume of NaOH solution required (ml)	Mean volume of NaOH solution (ml)	Strength of HCl solution
1	25	19.5	19.5	0.8823 (N/20)
2	25	19.5		

### Calculation:

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$\begin{aligned} \text{Strength of HCl solution (S}_1) &= \frac{19.5 \times 1.1312 \text{ (N/20)}}{25} \\ &= 0.8823 \text{ (N/20)} \end{aligned}$$

$V_1$  = Vol. of HCl solution

$V_2$  = Vol. of NaOH solution

$S_2$  = Strength of NaOH solution

**Table-3: Titration using phenolphthalein as indicator**

No of titration	Volume of supplied mixture solution (ml)	Volume of HCl solution required (ml)	Mean volume of HCl solution (ml)
1	25	9.5	9.5
2	25	9.4	
3	25	9.6	

**Table-4: Titration using methyl orange as indicator**

No of titration	Volume of supplied mixture solution (ml)	Volume of HCl solution required (ml)	Mean volume of HCl solution (ml)
1	25	23.5	23.5
2	25	23.4	
3	25	23.6	

### Calculation

#### Estimation of $\text{Na}_2\text{CO}_3$ solution:

Volume of HCl required for 1/2 mole  $\text{Na}_2\text{CO}_3$  solution = 9.5 ml

Volume of HCl required for 1mole  $\text{Na}_2\text{CO}_3$  solution =  $2 \times 9.5 \text{ ml} = 19 \text{ ml}$

$$\begin{aligned} 1000 \text{ ml } 1 \text{ (N) HCl solution} &\equiv 53 \text{ g of } \text{Na}_2\text{CO}_3 \\ 19 \text{ ml } 0.8823 \text{ (N/20)} &= \frac{53 \times 19 \times 0.8823}{1000 \times 20} \\ &= 0.0444 \text{ g of } \text{Na}_2\text{CO}_3 \end{aligned}$$

$$\begin{aligned} 25 \text{ ml supplied mixture sample solution} &\text{ contains } 0.0444 \text{ g of } \text{Na}_2\text{CO}_3 \\ 1000 &= 0.0444 \times 40 \\ &= 1.7769 \text{ g of } \text{Na}_2\text{CO}_3 \end{aligned}$$

$\therefore$  The amount of  $\text{Na}_2\text{CO}_3$  in the supplied mixture sample = 1.7769 g /lt

#### Estimation of $\text{NaHCO}_3$ :

Volume of HCl required for  $\text{NaHCO}_3$  solution present in the mixture =  $23.5 - 2 \times 9.5 = 4.5 \text{ ml}$

1000 ml 1(N) HCl solution  $\equiv 84 \text{ g of } \text{NaHCO}_3$

$$4.5 \text{ ml } 0.8823 \text{ (N/20)} \equiv \frac{84 \times 4.5 \times 0.8823}{1000 \times 20} = 0.01667 \text{ g of } \text{NaHCO}_3$$

Prepared by A.B Sir, Surendranath College, Please contact for any clarification.

25 ml supplied solution contains 0.01667 g  $\text{NaHCO}_3$

$$\begin{aligned} 1000 \text{ ml} \quad & \quad \quad \quad " \quad = \frac{0.01667 \times 1000}{25} \text{ g NaHCO}_3 \\ & \quad \quad \quad \quad \quad \quad \quad = 0.6668 \text{ g NaHCO}_3 \end{aligned}$$

$\therefore$  Therefore,  $\text{NaHCO}_3$  present in supplied solution = 0.6668 g/L

*Dr. A. B. Biswas*

**Dr. Apurba Biswas**  
**Assistant Professor of Chemistry**  
**Surentranath College**  
**Ph.D, Research Associate (C.U)**  
**Postdoc (South Korea & China)**

## Estimation of total hardness of water sample by complexometric EDTA titration

Date:

**Principle:**

Hardness of water is of two types- temporary hardness and permanent hardness. Temporary hardness is due to dissolved bicarbonates and permanent hardness is due to dissolved chlorides and sulphates, mainly of  $\text{Ca}^{2+}$ ,  $\text{Mg}^{2+}$  and  $\text{Fe}^{2+}$ . Hardness is expressed in terms of ppm of  $\text{CaCO}_3$  (i.e. mg of  $\text{CaCO}_3$  per 1000 ml of sample water).

Total hardness may be estimated complexometrically by titrating a known volume of the water sample with a standard EDTA solution in  $\text{NH}_4\text{Cl-NH}_3$  buffer medium ( $\text{pH}$  10) using Erichrome Black- T (EBT) indicator.

EDTA, ethylenediamine tetraacetic acid, is commercially obtained as its disodium salt,  $\text{Na}_2\text{H}_2\text{EDTA}\cdot 2\text{H}_2\text{O}$ . If it reacts with metal ions ( $\text{M}^{2+} = \text{Ca}^{2+}, \text{Mg}^{2+}, \text{Fe}^{2+}$ ) to form stable 1:1 complexes.



$$\therefore \text{EDTA} \equiv \text{M}^{2+} \equiv \text{CaCO}_3$$

$$\therefore 1000 \text{ ml } 1 \text{ (M) EDTA} \equiv 100 \text{ g of CaCO}_3$$

**Experimental Data**

**Table-1:** Standardization of supplied EDTA solution with the supplied standard (M/100) zinc acetate solution

No of titration	Volume of zinc acetate solution (ml)	Volume of EDTA solution required (ml)	Mean volume of EDTA solution (ml)	Strength of EDTA solution
1	25	27.0	27	0.9259 (M/100)
2	25	26.9		
3	25	27.1		

**Calculation:**

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$\begin{aligned} \text{Strength of EDTA solution (S}_2) &= \frac{V_1 \times S_1}{V_2} \\ &= \frac{25 \times \text{M}/100}{27} \\ &= 0.9259 \text{ (M/100)} \end{aligned}$$

$V_1$  = Vol. of zinc acetate solution

$S_1$  = Strength of zinc acetate solution

$V_2$  = Vol. of EDTA solution

**Table-2:** Estimation Hardness of water sample

No of titration	Volume of water sample (ml)	Volume of EDTA solution required (ml)	Mean volume of EDTA solution (ml)
1	50	7	7
2	50	7	

**Calculation:**

$$\begin{aligned} 1000 \text{ ml } 1 \text{ (M) EDTA solution} &\equiv 100 \text{ g of CaCO}_3 \\ 7 \text{ " } 0.9259 \text{ (M/100)} &= \frac{100 \times 7 \times 0.9259}{1000 \times 100} \text{ g of CaCO}_3 \\ &= 6.4813 \times 10^{-3} \text{ g of CaCO}_3 \end{aligned}$$

$$\begin{aligned} 50 \text{ ml water sample contains } &6.4813 \times 10^{-3} \text{ g of CaCO}_3 \\ 10^6 &= \frac{6.4813 \times 10^{-3} \times 10^6}{50} \end{aligned}$$

Total hardness of supplied water = 129.626 ppm

*Dr. Apurba Biswas*

**Dr. Apurba Biswas**  
Assistant Professor of Chemistry  
Surendranath College  
Ph.D, Research Associate (C.U)  
Postdoc (South Korea & China)