

**Estimation of sodium carbonate and sodium hydrogen carbonate present in a mixture**

Date:

**Principle:**

Neutralisation of Na<sub>2</sub>CO<sub>3</sub> solution by strong acid (HCl) occurs in two steps:



The acid base indicator is to be so selected that its P<sup>H</sup> range for the colour change coincides with the sudden sharp change of P<sup>H</sup> at the equivalence point. So at the first neutralisation point ( PH = 8.3), phenolphthalein (P<sup>H</sup> range 3.1-10) shows its colour change from pink to colourless. At this stage Na<sub>2</sub>CO<sub>3</sub> consumes only half the amount of HCl required for complete neutralisation. If methyl orange ( P<sup>H</sup> range = 3.1 to 4.4) is added to this titrated solution and the titration with HCl is continued upto second equivalence point, then this titre value corresponds to the amount of HCl required to convert NaHCO<sub>3</sub> to NaCl (i.e. NaHCO<sub>3</sub> derived from Na<sub>2</sub>CO<sub>3</sub> plus the amount of NaHCO<sub>3</sub> present in the original mixture).

1000ml 1(N) HCl  $\equiv$  53 g of Na<sub>2</sub>CO<sub>3</sub>  $\equiv$  84 g of NaHCO<sub>3</sub>

**Experimental Data**

**Table-1:** Standardization of supplied NaOH solution with the supplied standard (N/20) oxalic acid solution

No of titration	Volume of oxalic acid solution (ml)	Volume of NaOH solution required (ml)	Mean volume of NaOH solution (ml)	Strength of NaOH solution
1	25	22.1	22.1	1.1312 (N/20)
2	25	22.1		

**Calculation:**

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$\text{Strength of NaOH solution (S}_2) = \frac{25 \times \text{N}/20}{22.1} = 1.1312 \text{ (N/20)}$$

Dr. Apurba Biswas  
 Assistant Professor of Chemistry  
 Surendranath College  
 Ph.D, Research Associate (C.U)  
 Postdoc (South Korea & China)

$V_1$  = Vol. of oxalic acid solution  
 $S_1$  = Strength of oxalic acid solution  
 $V_2$  = Vol. of NaOH solution

**Table-2:** Standardization of supplied HCl solution with the standardized (N/20) NaOH solution

No of titration	Volume of HCl solution (ml)	Volume of NaOH solution required (ml)	Mean volume of NaOH solution (ml)	Strength of HCl solution
1	25	19.5	19.5	0.8823 (N/20)
2	25	19.5		

Write the exact reading that you got in your experiment and calculate accordingly

**Calculation:** Applying  $V_1 \times S_1 = V_2 \times S_2$

$$\begin{aligned} \text{Strength of HCl solution (S}_1) &= \frac{19.5 \times 1.1312 \text{ (N/20)}}{25} \\ &= 0.8823 \text{ (N/20)} \end{aligned}$$

$V_1$  = Vol. of HCl solution

$V_2$  = Vol. of NaOH solution

$S_2$  = Strength of NaOH solution

**Table-3: Titration using phenolphthalein as indicator**

No of titration	Volume of supplied mixture solution (ml)	Volume of HCl solution required (ml)	Mean volume of HCl solution (ml)
1	25	9.5	9.5
2	25	9.4	
3	25	9.6	

**Table-4: Titration using methyl orange as indicator**

No of titration	Volume of supplied mixture solution (ml)	Volume of HCl solution required (ml)	Mean volume of HCl solution (ml)
1	25	23.5	23.5
2	25	23.4	
3	25	23.6	

**Calculation**

**Estimation of  $\text{Na}_2\text{CO}_3$  solution:**

Volume of HCl required for 1/2 mole  $\text{Na}_2\text{CO}_3$  solution = 9.5 ml

Volume of HCl required for 1 mole  $\text{Na}_2\text{CO}_3$  solution =  $2 \times 9.5 \text{ ml} = 19 \text{ ml}$

$$\begin{aligned} 1000 \text{ ml } 1 \text{ (N) HCl solution} &\equiv 53 \text{ g of } \text{Na}_2\text{CO}_3 \\ 19 \text{ ml } 0.8823 \text{ (N/20)} &= \frac{53 \times 19 \times 0.8823}{1000 \times 20} \\ &= 0.0444 \text{ g of } \text{Na}_2\text{CO}_3 \end{aligned}$$

$$\begin{aligned} 25 \text{ ml supplied mixture sample solution} &\text{ contains } 0.0444 \text{ g of } \text{Na}_2\text{CO}_3 \\ 1000 &= 0.0444 \times 40 \\ &= 1.7769 \text{ g of } \text{Na}_2\text{CO}_3 \end{aligned}$$

$\therefore$  The amount of  $\text{Na}_2\text{CO}_3$  in the supplied mixture sample = 1.7769g/lt.

**Estimation of  $\text{NaHCO}_3$ :**

Volume of HCl required for  $\text{NaHCO}_3$  solution present in the mixture =  $23.5 - 2 \times 9.5 = 4.5 \text{ ml}$

1000 ml 1(N) HCl solution  $\equiv 84 \text{ g of } \text{NaHCO}_3$

$$4.5 \text{ ml } 0.8823 \text{ (N/20)} \equiv \frac{84 \times 4.5 \times 0.8823}{1000 \times 20} = 0.01667 \text{ g of } \text{NaHCO}_3$$

25 ml supplied solution contains 0.01667g  $\text{NaHCO}_3$

$$1000 \text{ ml} \equiv \frac{0.01667 \times 1000}{25} \text{ g } \text{NaHCO}_3$$

$$= 0.6668 \text{ g } \text{NaHCO}_3$$

$\therefore$  Therefore,  $\text{NaHCO}_3$  present in supplied solution = 0.6668g/L

Write the exact reading that you got in your experiment and calculate accordingly



## Experiment 2

### Estimation of oxalic acid by titrating it with $\text{KMnO}_4$

Date:

#### Principle:

Oxalic acid is estimated by titrating with standardised  $\text{KMnO}_4$  solution in 2(N)  $\text{H}_2\text{SO}_4$  medium at 70-80  $^\circ\text{C}$ . In dilute  $\text{H}_2\text{SO}_4$  acid medium  $\text{MnO}_4^-$  quantitatively oxidises  $\text{C}_2\text{O}_4^{2-}$  to  $\text{CO}_2$  and itself reduced to  $\text{Mn}^{2+}$ ,



$\text{KMnO}_4$  solution may be standardised against standard oxalic acid solution in 2(N)  $\text{H}_2\text{SO}_4$  medium at 70-80  $^\circ\text{C}$ . Purple coloured  $\text{KMnO}_4$  serves as self indicator. Its strengths may be calculated using the relation:

$$V_{\text{MnO}_4^-} \times S_{\text{MnO}_4^-} = V_{\text{C}_2\text{O}_4^{2-}} \times S_{\text{C}_2\text{O}_4^{2-}}$$

#### Experimental Data

(a) Table 1: Standardization of supplied  $\sim(\text{N}/20)$   $\text{KMnO}_4$  solution

No of titration	Volume of standard oxalic acid taken (ml)	Volume of $\text{KMnO}_4$ solution required (ml)	Mean volume of $\text{KMnO}_4$ solution required (ml)	Strength of $\text{KMnO}_4$ solution
1	25	25.6	25.6	0.9858(N/20)
2	25	25.7		
3	25	25.5		

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$25 \times 1.0095 (\text{N}/20) = 25.6 \times S_2$$
$$S_2 (\text{Strength of } \text{KMnO}_4) = \frac{25 \times 1.0095}{25.6} (\text{N}/20)$$
$$= 0.9858 (\text{N}/20)$$

$V_1$  = Vol. of oxalic acid

$S_1$  = Strength of oxalic acid

$V_2$  = Vol. of  $\text{KMnO}_4$

(b) Table-2: Estimation of oxalic acid

No of titration	Volume of supplied solution (ml)	Volume of $\text{KMnO}_4$ solution required (ml)	Mean volume of $\text{KMnO}_4$ solution required (ml)
1	25	12.3	12.3
2	25	12.4	
3	25	12.2	

#### Calculation:

1000 ml 1(N)  $\text{KMnO}_4$  solution  $\equiv$  63.03 g of oxalic acid

$$12.3 \text{ ml } 0.9858 (\text{N}/20) \text{ " } = \frac{63.03 \times 12.3 \times 0.9858}{1000 \times 20} = 0.0382 \text{ g of oxalic acid}$$

Dr. Apurba Biswas  
Assistant Professor of Chemistry  
Surenranath College  
Ph.D, Research Associate (C.U)  
Postdoc (South Korea & China)

**Write the exact reading that you got in your experiment and calculate accordingly**

25 ml supplied solution contains 0.0382 g of oxalic acid

$$\begin{aligned} 1000 \text{ ml} & \quad " & = \frac{0.0382 \times 1000}{25} \text{ g of oxalic acid} \\ & & = 1.528 \text{ g of oxalic acid} \end{aligned}$$

∴ Therefore, oxalic acid present in supplied solution = 1.528 g/L

*Dr. Apurba Biswas*

Dr. Apurba Biswas  
Assistant Professor of Chemistry  
Surendranath College  
Ph.D, Research Associate (C.U)  
Postdoc (South Korea & China)


## Experiment 3

**Estimation of Fe (II) ions by titrating it with K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> using internal indicator.**

Date:

**Principle:**

In acid medium Fe<sup>II</sup> in a solution may be estimated by direct titration with a standard solution of K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> in presence of either phosphoric acid (H<sub>3</sub>PO<sub>4</sub>) or ammonium bifluoride (NH<sub>4</sub>HF<sub>2</sub>) using barium diphenylamine sulphonate (BDS) as indicator. Under this condition K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> quantitatively oxidises Fe<sup>2+</sup> to Fe<sup>3+</sup>:



$$[\text{Cr}_2\text{O}_7^{2-}]/6 \equiv [\text{Fe}^{3+}]$$

$$1000\text{ml } 1(\text{N}) \text{ K}_2\text{Cr}_2\text{O}_7 \text{ solution} \equiv 55.847 \text{ g of Fe}^{2+}$$

H<sub>3</sub>PO<sub>4</sub> (or NH<sub>4</sub>HF<sub>2</sub>) stabilises Fe<sup>3+</sup> by complex formation, which is essential for indicator action of BDS.

**Experimental Data:****Determination of Fe<sup>2+</sup>:**

No. of titrations	Vol. of supplied Mohr's salt solution taken (ml)	Burette reading of standard K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> solution (ml)			
		Initial	Final	Difference	Mean
1	25	0	25.2	25.2	25.2
2	25	0	25.1	25.1	
3	25	0	25.3	25.3	

**Calculation:**

$$1000\text{ml } 1(\text{N}) \text{ K}_2\text{Cr}_2\text{O}_7 \text{ solution} \equiv 55.847 \text{ g of Fe}^{2+}$$

$$25.2\text{ml } 0.9915 \text{ (N/20) " } \equiv \frac{55.847 \times 25.2 \times 0.9915}{1000 \times 20} = 0.069977 \text{ g of Fe}^{2+}$$

$$25 \text{ ml supplied solution contains } 0.069977 \text{ g Fe}^{2+}$$

$$1000 \text{ ml " } = \frac{0.069977 \times 1000}{25} \text{ g Fe}^{2+}$$

$$= 2.7988 \text{ g Fe}^{2+}$$

$$\therefore \text{Therefore, Fe}^{2+} \text{ present in supplied solution} = 2.7988 \text{ g/L}$$


**Experiment 4**

**Estimation of Fe(II) and Fe(III) in a given mixture using K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution**

Date:


**Principle:**

In acid medium, direct titration of the Fe<sup>II</sup> + Fe<sup>III</sup> mixture with standard K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution gives the titre value (V<sub>1</sub>) corresponding to the amount of Fe<sup>2+</sup>. To estimate Fe<sup>3+</sup> present in the mixture, it is first reduced to Fe<sup>2+</sup> with highly pure Al-foil in 6(N) HCl medium.



Acidity of the solution before titration is adjusted to 2(N) by diluting with distilled water. The resulting solution is then titrated with the same standard K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution. This titre value (V<sub>2</sub>) corresponds to the total iron [Fe<sup>2+</sup> + Fe<sup>3+</sup>]. The difference (V<sub>2</sub>-V<sub>1</sub>) corresponds to the amount of Fe<sup>3+</sup>.

In acid medium K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> quantitatively oxidises Fe<sup>2+</sup> to Fe<sup>3+</sup>:



**Experimental Data:**

**(i) Determination of Fe<sup>2+</sup>:**

No. of titrations	Vol. of supplied solution taken (ml)	Burette reading of standard K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> solution (ml)			
		Initial	Final	Difference	Mean
1	25	0	25.2	25.2	25.2
2	25	0	25.1	25.1	
3	25	0	25.3	25.3	

**(ii) Determination of total iron (Fe<sup>2+</sup> + Fe<sup>3+</sup>):**

No. of titrations	Volume of supplied solution (ml)	Burette reading of standard K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> solution (ml)			
		Initial	Final	Difference	Mean
1	25	0	49.2	49.2	49.2
2	25	0	49.1	49.1	
3	25	0	49.3	49.3	

**Calculation:**

1000ml 1(N) K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution  $\equiv$  55.847 g of Fe<sup>2+</sup>

25.2 ml 0.9915 (N/20) "  $\equiv \frac{55.847 \times 25.2 \times 0.9915}{1000 \times 20} = 0.069977 \text{ g of Fe}^{2+}$

Dr. Apurba Biswas  
Assistant Professor of Chemistry  
Surendranath College  
Ph.D, Research Associate (C.U)  
Postdoc (South Korea & China)

Write the exact reading that you got in your experiment and calculate accordingly

25 ml supplied solution contains 0.069977 g Fe<sup>2+</sup>

$$\begin{aligned} 1000 \text{ ml} & \quad " & = \frac{0.069977 \times 1000}{25} \text{ g Fe}^{2+} \\ & & = 2.7988 \text{ g Fe}^{2+} \end{aligned}$$

∴ Therefore, Fe<sup>2+</sup> present in supplied solution = 2.7988 g/L

The difference (49.2-25.2) = 24 ml 0.9915 K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution gives the amount of Fe<sup>3+</sup>.

1000ml 1(N) K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution ≡ 55.847 g of Fe<sup>3+</sup>

$$24 \text{ ml } 0.9915 \text{ (N/20) } " \quad \equiv \frac{55.847 \times 24 \times 0.9915}{1000 \times 20} = 0.06643 \text{ g of Fe}^{3+}$$

25 ml supplied solution contains 0.066643 g Fe<sup>3+</sup>

$$\begin{aligned} 1000 \text{ ml} & \quad " & = \frac{0.06643 \times 1000}{25} \text{ g Fe}^{3+} \\ & & = 2.6572 \text{ g Fe}^{3+} \end{aligned}$$

∴ Therefore, Fe<sup>3+</sup> present in supplied solution = 2.6572 g/L

Dr. Apurba Biswas  
Assistant Professor of Chemistry  
Surendranath College  
Ph.D, Research Associate (C.U)  
Postdoc (South Korea & China)





Experiment 5

**Estimation of Cu (II) ions iodometrically using Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>**

**Date:**

**Principle:**

Copper (II) is iodometrically estimated by treating the solution containing Cu<sup>2+</sup> with an excess of KI solution, when sparingly soluble cuprous iodide, (Cu<sub>2</sub>I<sub>2</sub>) is precipitated with liberation of equivalent amount of iodine. Cu(II) is estimated by titrating the liberated iodine with a standardised solution of sodium thiosulphate.



$$1000\text{ml } 1(\text{N}) \text{ thiosulfate solution} \equiv 63.546 \text{ g of Cu}^{2+}$$

**Experimental Data:**

**(a) Standardization of supplied ~ (N/20) thiosulfate solution:**

No. of titrations	Vol. of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub> solution taken (ml)	Burette reading of thiosulphate solution (ml)			
		Initial	Final	Difference	Mean
1	25	0	26.4	26.4	26.3
2	25	0	26.3	26.3	
3	25	0	26.4	26.2	

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$25 \times 0.9915 (\text{N}/20) = 26.3 \times S_2$$

$$S_2 (\text{Strength of thiosulphate solution}) = \frac{25 \times 0.9915}{26.3} (\text{N}/20)$$

$$= 0.9425 (\text{N}/20)$$

$V_1$  = Vol. of K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution

$S_1$  = Strength of K<sub>2</sub>Cr<sub>2</sub>O<sub>7</sub> solution

$V_2$  = Vol. of thiosulfate solution

**(b) Estimation of Cu(II) by standardized thiosulfate solution:**

No. of titrations	Volume of Cu <sup>2+</sup> solution taken (ml)	Burette reading of thiosulphate solution (ml)			
		Initial	Final	Difference	Mean
1	25	0	30.1	30.1	30.1
2	25	0	30.1	30.1	
3	25	0	30.1	30.1	

Dr. Apurba Biswas  
 Assistant Professor of Chemistry  
 Surendranath College  
 Ph.D, Research Associate (C.U)  
 Postdoc (South Korea & China)

**Write the exact reading that you got in your experiment and calculate accordingly**

### Calculation

1000 ml 1(N) thiosulphate solution  $\equiv$  63.546 g of  $\text{Cu}^{2+}$

$$30.1 \text{ ml } 0.9425 \text{ (N/20)} \quad \equiv \quad \frac{63.546 \times 30.1 \times 0.9425}{1000 \times 20} = 0.0901 \text{ g } \text{Cu}^{2+}$$

25 ml stock solution contains 0.0901 g  $\text{Cu}^{2+}$

$$1000 \text{ ml} \quad = \frac{0.0901 \times 1000}{25} = 3.604 \text{ g } \text{Cu}^{2+}$$

$\therefore$  Therefore,  $\text{Cu}^{2+}$  present in supplied solution = 3.604 g/L

*Dr. Apurba Biswas*

Dr. Apurba Biswas  
Assistant Professor of Chemistry  
Surendranath College  
Ph.D, Research Associate (C.U)  
Postdoc (South Korea & China)

Write the exact reading that you got in your experiment and calculate accordingly

## Experiment 6

### Estimation of water of crystallization in Mohr's salt by titrating with $\text{KMnO}_4$

Date:

#### Principle:

Normality of ferrous ammonium sulphate can be determined by directly titrating it against standard (N/20)  $\text{KMnO}_4$  solution.

Equivalent weight = Strength/ Normality

$$\therefore \text{Equivalent weight of Mohr's salt} = \frac{\text{Strength of Mohr's salt (g/l)}}{\text{Normality of Mohr's salt}}$$

Substituting the value of strength and value of normality as calculated above, the equivalent weight of Mohr's salt can be calculated. Suppose it comes out to be E.


Since equivalent weight of Mohr's salt is equal to its molecular mass, therefore, molecular mass of Mohr's salt also equal to E.

Theoretical molecular mass of Mohr's salt,  $\text{FeSO}_4 \cdot (\text{NH}_4)_2\text{SO}_4 \cdot n\text{H}_2\text{O} = 284 + 18n$

$$284 + 18n = E \text{ where } n \text{ can be calculated.}$$

$$\therefore \text{No. of molecules of water of crystallisation, } n = \frac{E - 284}{18}$$

Ionic equation:



#### Experimental Data:

##### (a) Standardization of supplied ~ (N/20) $\text{KMnO}_4$ solution:

No of titration	Volume of standard oxalic acid taken (ml)	Volume of $\text{KMnO}_4$ solution required (ml)	Mean volume of $\text{KMnO}_4$ solution required (ml)	Strength of $\text{KMnO}_4$ solution
1	25	23.2	23.2	1.0830 (N/20)
2	25	23.1		
3	25	23.3		

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$\begin{aligned} 25 \times 1.005 \text{ (N/20)} &= 23.2 \times S_2 \\ S_2 \text{ (Strength of } \text{KMnO}_4) &= \frac{25 \times 1.005}{23.2} \text{ (N/20)} \\ &= 1.0830 \text{ (N/20)} \end{aligned}$$

$$\begin{aligned} V_1 &= \text{Vol. of oxalic acid} \\ S_1 &= \text{Strength of oxalic acid} \\ V_2 &= \text{Vol. of } \text{KMnO}_4 \end{aligned}$$

Write the exact reading that you got in your experiment and calculate accordingly

**(b) Standardization of Mohr's salt solution by standardized  $\text{KMnO}_4$  solution:**

No. of titrations	Vol. of Mohr's salt solution taken (ml)	Burette reading of standardized $\text{KMnO}_4$ solution (ml)			
		Initial	Final	Difference	Mean
1	25	0	22.6	22.6	22.6
2	25	0	22.6	22.6	
3	25	0	22.6	22.6	

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$25 \times S_1 = 22.6 \times 1.0830 \text{ (N/20)}$$

$$\begin{aligned} S_1 \text{ (Strength of Mohr's salt solution)} &= \frac{22.6 \times 1.0830}{25} \text{ (N/20)} \\ &= 0.9790 \text{ (N/20)} \\ &= 0.04895 \text{ (N)} \end{aligned}$$

$V_1$  = Vol. of Mohr's salt solution

$V_2$  = Vol. of  $\text{KMnO}_4$  solution

$S_2$  = Strength of  $\text{KMnO}_4$  solution)

**Calculation of equivalent weight of Mohr's salt:**

$$\begin{aligned} \text{Equivalent weight of Mohr's salt, } E &= \frac{\text{Strength of Mohr's salt (g/l)}}{\text{Normality of Mohr's salt}} \\ &= \frac{19.6}{0.04895} \end{aligned}$$

**No. of molecules of water of crystallisation, n:**

$$\therefore 284 + 18n = \frac{19.6}{0.04895}$$

$$18n = 400.40858 - 284$$

$$n = 6.467$$


**Experiment 7 (Not in the syllabus)**

**Estimation of Fe(II) using standardized KMnO<sub>4</sub> solution**

Date:

**Principle:**

In dilute H<sub>2</sub>SO<sub>4</sub> acid medium, KMnO<sub>4</sub> quantitatively oxidises Fe<sup>2+</sup> to Fe<sup>3+</sup>:



$$\therefore \text{MnO}_4^- \equiv 5\text{Fe}^{2+}$$

$$\therefore \text{MnO}_4^-/5 \equiv \text{Fe}^{2+}$$

Mohr's salt (Fe<sup>2+</sup>) solution may be estimated by titrating it, in 2(N) H<sub>2</sub>SO<sub>4</sub> medium with a standardised solution of KMnO<sub>4</sub> at room temperature in presence of H<sub>3</sub>PO<sub>4</sub>.

$$1000 \text{ ml } 1(\text{N}) \text{ KMnO}_4 \equiv 55.847 \text{ g of Fe}^{2+}$$

**Experimental Data:**

**(a) Standardization of supplied ~ (N/20) KMnO<sub>4</sub> solution:**

No of titration	Volume of standard oxalic acid taken (ml)	Volume of KMnO <sub>4</sub> solution required (ml)	Mean volume of KMnO <sub>4</sub> solution required (ml)	Strength of KMnO <sub>4</sub> solution
1	25	23.2	23.2	1.0830 (N/20)
2	25	23.1		
3	25	23.3		

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$25 \times 1.005 (\text{N}/20) = 23.2 \times S_2$$

$$S_2 (\text{Strength of KMnO}_4) = \frac{25 \times 1.005}{23.2} (\text{N}/20) = 1.0830 (\text{N}/20)$$

V<sub>1</sub> = Vol. of oxalic acid

S<sub>1</sub> = Strength of oxalic acid

V<sub>2</sub> = Vol. of KMnO<sub>4</sub>

**(b) Standardization of Mohr's salt solution/ Estimation of Fe(II) in Mohr's salt by standardized KMnO<sub>4</sub> solution:**

No. of titrations	Vol. of Mohr's salt solution taken (ml)	Burette reading of standardized KMnO <sub>4</sub> solution (ml)			
		Initial	Final	Difference	Mean
1	25	0	22.6	22.6	22.6
2	25	0	22.6	22.6	
3	25	0	22.6	22.6	

Applying  $V_1 \times S_1 = V_2 \times S_2$

$$25 \times S_1 = 22.6 \times 1.0830 (\text{N}/20)$$

$$S_1 (\text{Strength of Mohr's salt solution}) = \frac{22.6 \times 1.0830}{25} (\text{N}/20) = 0.9790 (\text{N}/20)$$

V<sub>1</sub> = Vol. of Mohr's salt solution

V<sub>2</sub> = Vol. of KMnO<sub>4</sub> solution

S<sub>2</sub> = Strength of KMnO<sub>4</sub> solution

Write the exact reading that you got in your experiment and calculate accordingly

### Calculation for the estimation of Fe(II) in Mohr's salt solution

1000 ml 1(N)  $\text{KMnO}_4$   $\equiv$  55.847 g of  $\text{Fe}^{2+}$

22.6 ml 1.0830 (N/20) "  $\equiv \frac{55.847 \times 22.6 \times 1.0830}{1000 \times 20} = 0.0683$  g of  $\text{Fe}^{2+}$

25 ml supplied Mohr's salt solution contains 0.0683 g  $\text{Fe}^{2+}$

1000 ml " "  $= \frac{0.0683 \times 1000}{25}$  g  $\text{Fe}^{2+} = 2.732$  g  $\text{Fe}^{2+}$

$\therefore$  Therefore,  $\text{Fe}^{2+}$  present in supplied Mohr's salt solution = 2.732 g/L

Dr. Apurba Biswas  
Assistant Professor of Chemistry  
Surendranath College  
Ph.D, Research Associate (C.U)  
Postdoc (South Korea & China)